

November 18, 2022

Zyno Medical has been notified of the use of IV administration sets for the Z-800 IV Infusion Pump that are not manufactured or approved for use with the pump by Zyno Medical. The Zyno IV administration sets are developed for specific use with the Z-800 Infusion pump. The devices are designed and manufactured with the corresponding IV administration sets for the device to ensure the flow rate accuracy as stated in the specifications.

The use of any other IV sets may cause improper device operation and result in inaccurate fluid delivery or other potential hazards. For reference, please note the following warning from the Z-800 Instructions for Use:

“Only use Zyno Medical provided IV sets with the Z-800 Infusion Pump. There are risks associated with using any IV Sets other than Zyno Medical IV Sets with this device. Zyno’s warranty for its device will be null and void and Zyno will assume no responsibility for any incidents that may occur if the device is not utilized strictly in accordance with its product labeling.”

Zyno Medical offers several IV administration set configurations that are designed specifically for use with the Z-800 IV Infusion Pump. Please contact your local Zyno Solutions representative or distribution partner for more information about the Zyno Medical IV administration sets. You can also contact us here: [support@zynomed.com](mailto:support@zynomed.com)

Sincerely,



Mark Coggin  
VP, Marketing  
Care Everywhere Group  
*Infutronic, Infutronic Solutions, Zyno Medical, Zyno Solutions, Care Everywhere*