



AGILIA® CONNECT VOLUMETRIC PUMP

MANUFACTURED BY FRESENIUS KABI

rightwaymed.com | 614.396.7721



Designed with usability engineering and smart technology in mind, and manufactured to be both durable and serviceable, the Agilia® Connect Volumetric Pump fits easily into every user's everyday workflow. It's your ally in safely and efficiently caring for patients in your healthcare facility.

Addressing important points in the infusion process, the Agilia Connect system of pumps, software, sets, and services, helps ensure everyone is armed with well-considered tools, aligned in their practice goals, and able to work as one big team.

The Agilia Connect Infusion System offering by Fresenius Kabi consists of four components:

- *Agilia Connect Pumps:* Volumetric Pump and Syringe Pump
- *Vigilant Software Suite:* a) Vigilant Master Med Dose Error Reduction Software, b) Vigilant Insight Continuous Quality Improvement Software, c) Vigilant Centerium Data Management Platform, d) Agilia Partner Service Software
- *Sets:* a) Volumat Line of Dedicated Sets, b) KabiAccess Non-Dedicated Disposables
- *Service:* A team of dedicated experts to help you get most out of your investment

SPECIFICATIONS

Dimensions: 135 x 190 x 170 mm

Weight: Approximately 2kg

Accuracy: ± 5%

Pumping Mechanism: 2nd generation linear peristaltic pumping system with software adjustments and compensations

Delivery Modes: mL/h modes: volume + flow rate, volume + time, flow rate + time, volume + time + rate, ramp-up / ramp down, sequential / intermittent, secondary / piggyback, drop/min

Battery: Li-ion Smart battery. WiFi disabled: > 8 h at 25 mL/h > 5 h at 1500 mL/h. WiFi enabled: 5 h at 25 mL/h, 4 h at 1500 mL/h.

Flow Rate: 0.1 - 1500 mL/h. 0.01 minimum increments instead of 0.1 mL/h is activated when the max patient weight is set to 20 kg in the drug library.

Volume: 0,1 - 9999 mL.

Prime: Manual or by using the purge/prime function of the pump

KVO: Automatic KVO rate of 1 mL/h (adjustable from 1 to 20 mL/h) when volume limit is reached

Scan For List
Of Fresenius Kabi
Products



ABOUT RIGHT WAY MEDICAL

Right Way Medical is a national concierge and consultant of medical products for alternate site, long-term care, and specialty pharmacy providers. We set the standard for transparency and responsiveness while delivering patient-ready infusion devices, disposables, biomedical services, and asset tracking software, Right Track.



Learn more at www.rightwaymed.com or 614.396.7721.