



CADD[®] SOLIS VIP

MANUFACTURED BY ICU MEDICAL

rightwaymed.com | 614.396.7721



The CADD[®]-Solis VIP (variable infusion profile) Ambulatory Infusion Pump is designed to facilitate patient care for a variety of adult and pediatric patients and clinical care areas. The pump can be programmed with a protocol configuration consisting of a therapy, qualifier and drug information.

The pump can deliver via the following delivery profiles: PCA, continuous, intermittent, variable stepped rate, and tapered infusions for analgesia, chemotherapy, hydration, antibiotics, immune globulin and TPN.

CADD[™]-Solis Medication Safety Software Administrator Program: Administrator software is a Windows[®] based application used by select staff members to create and manage the standardized pump protocols. The application is used to securely program the CADD[®]-Solis VIP pump saving Pharmacist and Nursing time in manually programming the infusion pump. Infusion data is collected for additional pump documentation, reporting and CQI data analysis.

SPECIFICATIONS

Dimensions: 4 in. x 5 in. x 1.6 in

Weight: 21 oz. including four AA batteries

Infusion Device: Multi-therapy ambulatory infusion pump

Pumping Mechanism: Linear peristaltic

Delivery Modes: Library programming (continuous, Intermittent, PCA, Taper, Step)

KVO Rate: Programmable, 0-10 mL/hr. (all therapies except PCA), PCA mode only: 0.1 mL/hr. – only when delivery limit reached

Accuracy: ±6%

Battery: 4 x AA batteries, rechargeable battery pack, AC power

Continuous Flow Rate: 0.1-250 mL/hr with standard set, 250-500 mL/hr. with high volume set, 0.1-100 mL/hr (PCA only)

Reservoir Volume Limits: 1.0 to 9,999 mL

Air Detector Alarm: LOW: single bubble > 2 mL, accumulated air: > 2 mL over 15 min. HIGH: Single bubble > 0.15 mL, accumulated air: > 1 mL over 15 min

Prime: 250 mL/hr standard set, 500 mL/hr high volume set

Preventative Maintenance Reminder: Programmable alert

Scan For List
Of ICU Medical
Products



ABOUT RIGHT WAY MEDICAL

Right Way Medical is a national concierge and consultant of medical products for alternate site, long-term care, and specialty pharmacy providers. We set the standard for transparency and responsiveness while delivering patient-ready infusion devices, disposables, biomedical services, and asset tracking software, Right Track.



Learn more at www.rightwaymed.com or 614.396.7721.