



SAPPHIRE MULTI-THERAPY

MANUFACTURED BY EITAN MEDICAL

rightwaymed.com | 614.396.7721



With a small, simplistic design, the Sapphire Multi-Therapy Infusion Pump is intuitive and easy to use. The large, color touch screen is reminiscent of a mobile application allowing for efficient programming, and is aligned with feature rich software for maximum accuracy & reliability.

INNOVATIVE

- Multi-therapy capability: TPN, PCA, continuous, multi-step, piggyback, intermittent, epidural, and PIEB
- Configurable settings and patient-specific parameters
- Easy to switch between pump and gravity infusion

SIMPLE

- Dual anti-free flow protection in pump and administration set
- Auto restart after downstream occlusion
- Resume option after shutdown during infusion
- On-screen alarm troubleshooting
- Rated for air and ground transport

ACCURATE

- Ability to change therapy without stopping device
- Large touchscreen providing ease of use during programming
- Four lock levels
- Highly durable with high protection against fluids, impact and power surges

SPECIFICATIONS

Dimensions: 5.63" H x 3.78" W x 1.93" D

Weight: 14.7 oz

Infusion Device: Volumetric, peristaltic

Pumping Mechanism: Single channel volumetric, with integral pressure sensor

Delivery Modes: Continuous (with and without secondary line), PCA, intermittent, multi-step, TPN, epidural (PCEA, epidural intermittent, PIEB)

KVO Rate: Up to 20 mL/h in increments of 0.1 mL/h

Accuracy: $\pm 2.5\%$ under normal conditions

Battery: Rechargeable Li-Ion battery, 24 hrs @125 mL/h with fully charged battery and backlight off

Flow Rate: 0.1-99.9 mL/h in increments of 0.1 mL/h, 100-999 mL/h in increments of 1 mL/h

Volume: 0.1 - 9999 mL in increments of 0.1 mL

Air-in-Line Sensor: Detects single and accumulated bubbles of selectable sizes varying from .02 - 0.5 mL

Prime: Manual or automatic priming

Scan For List
Of Eitan Medical
Products



ABOUT RIGHT WAY MEDICAL

Right Way Medical is a national concierge and consultant of medical products for alternate site, long-term care, and specialty pharmacy providers. We set the standard for transparency and responsiveness while delivering patient-ready infusion devices, disposables, biomedical services, and asset tracking software, Right Track.



Learn more at www.rightwaymed.com or 614.396.7721.