



VISTA BASIC

MANUFACTURED BY B. BRAUN

rightwaymed.com | 614.396.7721



The Vista basic was designed to ensure simplicity in all aspects of its use: from set installation, to programming, to special mode use and alarm correction. Based on a design that has sold over 100,000 units worldwide, the Vista basic provides high reliability packaged in one of the easiest pumps to use.

As a workhorse in the outpatient market, it offers a rugged and reliable design, all relevant therapies such as continuous, piggyback and program mode, and a selection of cost-effective, straight line gravity sets.

- Robust design of mechanical components ensures maximum durability
- Battery life of approximately 3 years before replacement is required
- Sets not made with DEHP to protect patients. Needle-free sets protect both the patient and the clinician

SPECIFICATIONS

Dimensions: 5.5" W x 9.4" H x 7.8" D

Weight: 6.8 lbs.

Pumping Mechanism: Linear peristaltic volumetric infusion pump

Delivery Modes: Continuous and Piggyback Modes, Program Mode for variable infusions, Ramp/Taper Mode for TPN infusions and Delay A Mode for postponing infusion start and finish

KVO Rate: 1 mL/hr. or selected rate (whichever is lower)

Accuracy: ± 5% measurement time 8hr, temp. 71.6° F

Battery: 7.2V, 1.2 Ah NiCd (rechargeable)

Battery Operating Time: 3.5 hours at 125 mL/hr.

Rate Range: 0.1 to 800.0 mL/hr. in 0.1 mL increments with Vista™ tubing

Volume to be Infused Range: 0.1 to 9999.9 mL.

Air-in-Line Alarm: Air bubble size 0.3mL* or 1.5mL/hr.** (can be change in service program- *0.1mL-0.3mL, *0.5 to 3.5 mL/hr)

Alarms: Air-in-Line, Low Battery, Plug in Now, Bag Empty, Downstream Occlusion, Door Open, Pump Stopped, KVO, Unit Malfunction

Scan For List
Of B. Braun
Products



ABOUT RIGHT WAY MEDICAL

Right Way Medical is a national concierge and consultant of medical products for alternate site, long-term care, and specialty pharmacy providers. We set the standard for transparency and responsiveness while delivering patient-ready infusion devices, disposables, biomedical services, and asset tracking software, Right Track.



Learn more at www.rightwaymed.com or 614.396.7721.