



Z-800F

MANUFACTURED BY INTUVIE

rightwaymed.com | 614.396.7721



The Z-800 series is an easy-to-use large-volume infusion system that is designed to support productive clinical workflow. Its intuitive programming and easy-to-load IV sets allow your team to start an infusion in a few steps. Built for long-term reliability, it's equipped with durable hardware that helps reduce pump downtime and repairs. With up to seven unique infusion modes, the Z-800 series is versatile for various infusion therapies and settings.

Anti-Free-Flow Protection

Prevent inadvertent fluid free flow when unloading IV set from the pump with pump-based and set-based anti-free-flow protection.

Access Protocol Mode

Eliminate repetitive programming for commonly used infusion regimens by saving protocols and retrieving them in the future to automatically program an infusion.

Nurse Notification Light

Observe, round, and assist patients easily with an external light connected to the pump and top of the IV pole—designed as add on feature & notification for pump alerts/alarms.

SPECIFICATIONS

Dimensions: 8.6" H x 5.7" W x 5.3" D

Weight: 7.0 lbs, 3.2 kg

Pumping Mechanism: Linear peristaltic

Infusion Modes: Continuous rate/volume, continuous time/volume, blood, ten step sequencing, TPN, access protocol, intermittent

Alarms & Warnings: Air-in-line, occlusion, door open, battery empty, no drip, no clamp, low battery, pump unattended, near end, infusion complete - KVO, due for service

Accuracy: ±5%

Battery: Rechargeable (9.6 V; 4.5 amp-hr.)

Battery Operating Time: 8 hours at 125 mL/h

Flow Rate: 1 - 1200 mL/h in 1 mL increment

Volume: 1 - 9999 mL in 1 mL increment

Air-in-Line Sensor: Ultrasonic, 14 level configurable

Priming Rate: 600 mL/h

Free Flow Protection: Pump based and set based in all models

Wireless Integration: No

Scan For List
Of Intuvie
Products



ABOUT RIGHT WAY MEDICAL

Right Way Medical is a national concierge and consultant of medical products for alternate site, long-term care, and specialty pharmacy providers. We set the standard for transparency and responsiveness while delivering patient-ready infusion devices, disposables, biomedical services, and asset tracking software, Right Track.



Learn more at www.rightwaymed.com or 614.396.7721.