

# Agilia<sup>®</sup> Connect

## Syringe Pump

Part of the **Agilia Connect** Infusion System



# Your high-acuity ally

Designed with usability engineering and smart technology in mind, and manufactured to be durable and serviceable, the Agilia® Connect Syringe Pump fits easily into every user's everyday workflow to accurately deliver critical care and high alert medications. From anesthesia to pediatrics, in the ER to the OR, it's your ally in promoting best practices in urgent moments.

## Clinically Focused on patient safety

In addition to supporting Vigilant Master Med, the Agilia Connect Syringe Pump has a set of built-in features and functions that help to keep safety front and center for your team.

## Technically Packed with powerful technology

The enhanced microprocessor and memory of the WiFi-enabled Agilia Connect Syringe Pump is able to handle regular uploads of Vigilant Master Med libraries and downloads of infusion data to Vigilant Insight, all without physically handling the pumps.

## Functionally Refined for intuitive operation

Human factors design has been applied to every aspect of the Agilia Connect Syringe Pump – from form to function – with a face, buttons, and programming that aligns to the Agilia Connect Volumetric Pump, to ease initial training and standardize everyday workflow across your entire pump fleet.

## Practically Engineered for reliable performance

The Agilia Connect Syringe Pump is built to perform with accuracy and – in conjunction with Agilia Partner – to maintain and repair with efficiency.

### Lightweight.

Less than

5 lbs

### Flexible.

Infusion rate from

0.1 to  
1200 mL/h

### Precise.

Infuses at

± 3%  
accuracy under  
standard conditions

### Reliable.

3 year  
Preventative  
Maintenance cycle



**Anti-bolus backflow pumping system** reduces volume after occlusion release

**Automatic detection** of syringe brands and sizes

**Protective syringe guard** prevents accidental dislodgement

**"Focus-friendly", incremental keypad** helps avoid key bounce and order of magnitude errors

**Visually-cued elements** like color-coded parts and pictogram button labels ease loading and speed set-up

**High-resolution screen** enables TALLman display of drug names and highlighting of concentrations, for at-a-glance visibility in all lighting

**Stackable design** occupies a smaller footprint making it easier to ambulate or transport patient

**Versatile pole clamp** attaches easily to an IV pole, bed rail, or even a table

**Wide range of delivery mode combinations**

- mL/h, mL/kg, body weight, mass volume over time, concentrations, bolus & loading dose, volume limit, KVO, standbys
- 2nd decimal point can be enabled or disabled

**Vigilant Master Med dose error reduction software**

- defines soft & hard limits for dose and flow rates
- accommodates a data set with up to 19 profiles that may include up to 3,800 different drug/therapies entries
- provides automatic weight-based dosing calculations

**Rugged ABS shell** provides shock protection

**High-performance battery**

- 3-year replacement cycle
- long-lasting between charges: > 6 hrs at 5mL/h with WiFi or > 11 hrs at 5mL/h without WiFi
- can recharge in < 6 hrs when turned off or < 20 hrs when on

**3-year preventative maintenance cycle** helps to minimize maintenance activities



### **Dedicated Experts**

Account Manager  
Project Manager  
Clinical Pharmacy Services Liaison  
Clinical Nurse Liaison  
IT Service Liaison

### **SCP certification, 10 years running**

Our technical team's commitment to annual certification by the global standard in field service quality – Service Capability & Performance (SCP) – whose audits cover a broad scope of practices, ensures installation and training, as well as ongoing service, are provided by individuals with deep knowledge.

### **Digital training tools**

E-learning resources bring Fresenius Kabi subject matter expertise to your teams anytime, anywhere, streamlining your new hire training and annual certifications.

## **Exceptionally Supported to drive best practices**

From initial implementation to everyday application, Fresenius Kabi clinical and technical experts work alongside your pharmacy, nursing, and biomed teams to put Agilia® Connect Syringe Pump to best use, getting everyone started quickly and ready to reach your infusion practice goals.

## **Globally Recognized as a proven partner**

Fresenius Kabi's medical device division is a leader throughout the world, with more than 35 years of infusion technology experience and an installed fleet of more than a million infusion pumps. Fresenius Kabi is also a steadily-growing force in the North American generic injectables space, having worked shoulder-to-shoulder with acute- and alternate-care institutions over the last decade to help ensure the availability and affordability of essential IV drugs for anesthesia and analgesia, critical care, and oncology.

SCP is a registered trademark of Service Strategies.



**FRESENIUS  
KABI**

caring for life

Fresenius Kabi, USA, LLC  
Three Corporate Drive  
Lake Zurich, Illinois 60047  
[www.fresenius-kabi.us/agiliconnect](http://www.fresenius-kabi.us/agiliconnect)

©2022 Fresenius Kabi, USA, LLC. All Rights Reserved.  
2575-ACIS-10-03/22