

WalkMed 350VL

WalkMed
Technologies

Ambulatory Infusion Pump



Ambulatory Infusion Pumps You Can Depend On

Safe, reliable ambulatory infusion equipment with a 30 year history of excellence.

Safety & Convenience

- 4 Step Programming
- Patient Lockout
- Volume Tracking
- Long Battery Life
- Small & Lightweight



WalkMed Reservoir Bags in 3 sizes:

- 65mL
- 150mL
- 250mL



WalkMed Admin Sets in 4 options:

- 40" Standard set
- 60" Standard set
- 0.2 Micron air eliminating filter set
- 0.2 Micron air eliminating filter set with back check valve

24 Hour Patient Helpline
WWW.WalkMed.net
Sales Email: Sales@WalkMed.net
Toll Free: (800) 578-0555

Precision • Quality • Performance

WalkMed 350VL

WalkMed
Technologies



Ambulatory Infusion Pump

24 Hour Patient Helpline:

- (800) 578-0555

Dimensions:

- 4.4 x 4.0 x 1.5 in (including front cover)

Weight:

- 12 oz (including battery and front cover)

Delivery Accuracy:

- +/- 5%

Infusion Mode:

- Continuous

Flow Rate:

- 0.10 - 19.99mL/hr (0.01mL Resolution)
- 20.0 - 30.0 mL/hr (0.10 mL Resolution)

Total Volume Limit:

- 1 - 1,999 mL

Reservoir Volume:

- 65, 150, 250 mL Bags

Power:

- 9 volt Battery (Standard)

Typical Battery Capacity:

- 1 mL/hr - 18.7 Days
- 10 mL/hr - 2.7 Days
- 30 mL/hr - 0.7 Days

Alarms/Alerts:

- System Idle • Downstream Occlusion • Low Battery
- Depleted Battery • Under-Delivery • Over-Delivery
- Clamp Bar Open • System Malfunction
- End of Infusion Notice • KVO/TVL Reached.



Precision • Quality • Performance